

# Product Information

## ELISA Tests for Grapevine Pathogens

Our ELISA reagents for the detection of grapevine viruses are provided in the convenient double antibody sandwich procedure (DAS-ELISA) for a working volume of 0.2 ml per well. They are offered either as «Complete kits», «Reagent sets» or «Individual reagents» in different sizes. For details, see [www.bioreba.com](http://www.bioreba.com) or our product catalogue. The following ELISA reagents are available:

### Group A:

Nepovirus

ArMV  
GFLV  
ArMV+GFLV  
RpRSV-ch  
RpRSV-g  
TBRV  
TRSV  
ToRSV  
ToRSV-Ch  
SLRSV

Arabis mosaic virus  
Grapevine fanleaf virus  
Arabis mosaic virus + Grapevine fanleaf virus  
Raspberry ringspot virus-ch  
Raspberry ringspot virus-g  
Tomato black ring virus  
Tobacco ringspot virus  
Tomato ringspot virus  
Tomato ringspot virus-Ch  
Strawberry latent ringspot virus

Sadwavirus (prev. Nepovirus)

### Group B:

Ampelo- and Closterovirus

GLRaV-1  
GLRaV-1+3  
GLRaV-2  
GLRaV-3  
GLRaV-4 strains  
GLRaV-4 strain 6  
GVA  
GFkV

Grapevine leafroll-assoc. virus 1  
Grapevine leafroll-assoc. virus 1+3  
Grapevine leafroll-assoc. virus 2  
Grapevine leafroll-assoc. virus 3  
Grapevine leafroll-assoc. virus generic 4 strains  
Grapevine leafroll-assoc. virus 4 strain 6  
Grapevine virus A  
Grapevine fleck virus

Vitivirus

Maculavirus

### Recommendations for sampling

The ELISA technique is an efficient method for the detection of these viruses in grapevine. However, the virus concentration varies considerably according to the tissue source, the meteorological conditions and thus, the time of the season. These facts have to be considered for obtaining reliable test results. For testing grapevine, a special extraction buffer «Grapevine» (Art. No. 110123) is used at a ratio of 1:10 (w/v).

The following recommendations for tissue sampling will help improve your ELISA test results:

#### For Group A (Nepo- and Sadwaviruses):

##### Preferable

Leaves from young shoots and juicy bark early in the growing season as well as bark (phloem) scrapings from mature canes during dormancy

##### Not recommended

Samples collected during the hot summer or old leaves in fall

#### For Group B (Ampelo-, Clostero-, Viti- and Maculaviruses):

##### Preferable

Well-developed mature leaves, especially petioles and veins, from the lower part of the plant late in the season as well as bark (phloem) scrapings from mature canes during dormancy

##### Not recommended

Young tissue early in the growing season

These recommendations fit best the conditions in the grape-growing areas in Switzerland and surroundings. Ideal sampling conditions may vary in climatic different growing conditions.

1/1

Version: 6 - 25.04.2017

Adaptations from last version: minor modifications.

#### Authorised distributor

##### In Australia:

For customer service, call 1300-735-292  
To fax an order, use 1800-067-639

To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

##### In New Zealand:

For customer service, call 0800-933-966  
To fax an order, use 0800-329-246

To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

**ThermoFisher**  
SCIENTIFIC

