



# The Fisherbrand

## Glassware collection

**ThermoFisher**  
SCIENTIFIC

# Volumetric Flasks

## Clear flask ground glass neck Class A

- ISO1042, DIN 12664 Pt I, BS1792 conform
- With cap
- Ring mark
- Borosilicate glass 3.3



Item code	Volume (mL)	Qty	Tolerance (mL)	Diameter (mm)	Height (mm)	Ground Neck	Price
FSBFB50141	5	5	± 0.025	22	69	10/19	\$25.00
FSBFB50143	10	5	± 0.025	28	89	10/19	\$21.00
FSBFB50145	20	5	± 0.040	36	111	10/19	\$23.00
FSBFB50147	25	5	± 0.040	38	110	10/19	\$23.00
FSBFB50149	50	5	± 0.060	49	139	12/21	\$27.00
FSBFB50151	100	5	± 0.100	61	170	12/22	\$31.00
FSBFB50153	200	2	± 0.150	76	210	14/23	\$22.00
FSBFB50155	250	2	± 0.150	81	220	14/23	\$25.00
FSBFB50157	500	2	± 0.250	95	260	19/26	\$29.00
FSBFB50159	1000	2	± 0.400	130	300	24/29	\$39.00
FSBFB50161	2000	1	± 0.600	147	370	29/32	\$29.00

## Amber flask ground glass neck Class A

- ISO1042, DIN 12664 Pt I, BS1792 conformity
- With plastic stopper
- Ring mark
- Borosilicate glass 3.3



Item code	Volume (mL)	Qty	Tolerance (mL)	Diameter (mm)	Height (mm)	Ground Neck	Price
FSBFB55533	10	5	± 0.025	28	89	10/19	\$219.00
FSBFB55535	20	5	± 0.040	36	111	10/19	\$99.00
FSBFB55537	25	5	± 0.040	38	110	10/19	\$119.00
FSBFB55539	50	5	± 0.060	49	139	12/21	\$129.00
FSBFB55541	100	5	± 0.100	61	170	12/21	\$159.00
FSBFB55543	200	2	± 0.150	76	210	14/23	\$99.00
FSBFB55545	250	2	± 0.150	81	220	14/23	\$99.00
FSBFB55547	500	2	± 0.250	95	260	19/26	\$119.00

# Measuring Cylinders

## Measuring cylinder clear Class A


- ISO4788, DIN12680-2, BS604 conformity
- Stable hexagonal foot
- With spout
- Delivered with conformity certificate
- Borosilicate glass 3.3



Item code	Volume (mL)	Qty	Tolerance (mL)	Diameter (mm)	Height (mm)	Price
FSBFB55670	10	2	± 0.10	15	135	\$19.00
FSBFB55671	25	2	± 0.25	21	160	\$16.00
FSBFB55672	50	2	± 0.50	26	195	\$17.00
FSBFB55673	100	2	± 0.50	31	240	\$18.00
FSBFB55674	250	2	± 1.00	40	320	\$24.00
FSBFB55675	500	2	± 2.50	53	380	\$35.00
FSBFB55676	1000	2	± 5.00	67	464	\$54.00

## Mixing cylinder ground neck Class B

- Borosilicate glass
- ISO 4788 conformity



Item code	Volume (mL)	Qty	Graduations (mL)	Tolerance (mL)	Diameter (mm)	Height (mm)	Ground Neck	Price
FSBFB55227	50	2	1	± 1.0	22	200	19/25	\$37.00
FSBFB55229	100	2	1	± 1.0	26	260	24/29	\$47.00
FSBFB55231	250	2	2	± 2.0	37	335	29/32	\$67.00
FSBFB55233	500	2	5	± 5.0	48	390	34/35	\$87.00
FSBFB55235	1000	2	10	± 10.0	61	470	45/40	\$127.00
FSBFB55237	2000	1	20	± 20.0	78	570	45/40	\$107.00

# Volmetric Pipettes

## Volumetric pipette Class AS

- ISO 648 DIN 12690 conform
- Total flow
- Borosilicate glass 3.3
- One mark

Item code	Volume (mL)	Qty	Tolerance (mL)	Length (mm)	Price
FSBFB50192	1	2	± 0.007	310	\$7.00
FSBFB50194	2	2	± 0.010	340	\$8.00
FSBFB50196	5	2	± 0.015	410	\$9.00
FSBFB50198	10	2	± 0.020	450	\$10.00
FSBFB50200	25	2	± 0.030	530	\$14.00
FSBFB50202	50	2	± 0.050	530	\$19.00




# Graduated Pipettes

## Graduated pipette type 1 Class AS

- ISO 835 conformity
- Soda lime glass
- Blue ceramic graduations

Item code	Volume (mL)	Qty	Graduations (mL)	Price
FSBFB50232	1	2	0.01	\$7.45
FSBFB50234	2	2	0.02	\$7.75
FSBFB50236	5	2	0.05	\$7.75
FSBFB50238	10	2	0.1	\$8.45
FSBFB50240	25	2	0.1	\$16.45



# Funnels

## Separating funnel conical form

- Manufactured from Borosilicate Glass
- ISO 4800 certified
- PTFE stopcock

Item code	Volume (mL)	Qty	Ground Neck	Price
FSBFB51741	100	1	19/26	\$49.00
FSBFB51742	250	1	24/29	\$59.00
FSBFB51743	500	1	24/29	\$69.00
FSBFB51744	1000	1	29/32	\$89.00



## Separating funnel pear form

- Manufactured from Borosilicate Glass
- ISO 4800 certified
- PTFE stopcock

Item code	Volume (mL)	Qty	Ground Neck	Price
FSBFB51750	50	1	19/26	\$49.00
FSBFB51751	100	1	19/26	\$55.00
FSBFB51752	250	1	24/29	\$69.00
FSBFB51753	500	1	24/29	\$75.00
FSBFB51754	1000	1	29/32	\$89.00



## Filter funnel sintered

- Funnels constructed of borosilicate glass.
- Funnel filters constructed of sintered glass
- Filters medium precipitates with grade 3 and grade 4 porosity

Item code	Volume (mL)	Porosity	Qty	Price
FSBFB68702	35	3	1	\$165.00
FSBFB68703	35	4	1	\$220.00
FSBFB68706	80	3	1	\$155.00
FSBFB68707	80	4	1	\$210.00
FSBFB68710	125	3	1	\$225.00
FSBFB68711	125	4	1	\$315.00
FSBFB68714	500	3	1	\$449.00



# Laboratory Bottles

- Permanent white enamel marking spot
- GL 45 screw thread
- Linerless
- Autoclavable (to 140°C)
- Borosilicate glass
- Polypropylene cap

Item code	Volume (mL)	Diameter (mm)	Height (mm)	Closure Size	Qty	Price
FSBFB33144	100	56	105	GL45	1	\$3.75
FSBFB33145	250	70	143	GL45	1	\$3.95
FSBFB33146	500	86	182	GL45	1	\$5.45
FSBFB33147	1000	101	230	GL45	1	\$7.75
FSBFB33148	2000	136	267	GL45	1	\$15.75
FSBFB33180	5000	181	330	GL45	1	\$75.00
FSBFB33181	10000	227	410	GL45	1	\$295.00



# Beakers

- Thick, heavy construction with beaded top
- Uniform wall thickness
- Durable, matte-finish marking area for use with an ordinary pencil
- Easy-to-read white graduated double scale on all sizes from 50 to 1000mL
- Borosilicate glass

## Low Form

Item code	Volume (mL)	Diameter (mm)	Height (mm)	Qty	Price
FSBFB10050	50	42	53	1	\$3.00
FSBFB100100	100	50	70	1	\$3.25
FSBFB100150	150	57	81	1	\$3.50
FSBFB100250	250	68	88	1	\$3.75
FSBFB100400	400	77	106	1	\$4.00
FSBFB100600	600	88	122	1	\$4.25
FSBFB1001000	1000	107	145	1	\$7.00



## Tall Form

Item code	Volume (mL)	Diameter (mm)	Height (mm)	Qty	Price
FSBFB102100	100	50	78	1	\$3.00
FSBFB102200	200	57	104	1	\$3.25
FSBFB102300	300	64	118	1	\$3.50
FSBFB102400	400	70	127	1	\$4.00
FSBFB102600	600	80	150	1	\$5.00
FSBFB1021000	1000	90	188	1	\$7.50




## Erlenmeyer Flasks

- Reinforced top tooled with a rounded finish
- Thick-walled body
- White ceramic enamel marking spot
- Long, tapered outside contour minimises chipping
- Borosilicate glass

### Narrow Neck

Item code	Volume (mL)	Diameter (mm)	Height (mm)	Qty	Price
FSBFB50025	25	39	65	1	\$3.50
FSBFB50050	50	50	78	1	\$3.75
FSBFB500125	125	66	108	1	\$3.00
FSBFB500250	250	82	130	1	\$3.50
FSBFB500500	500	102	174	1	\$4.00
FSBFB5001000	1000	128	213	1	\$7.00
FSBFB5002000	2000	160	275	1	\$13.00
FSBFB5004000	4000	208	355	1	\$35.00



### Wide Neck

Item code	Volume (mL)	Diameter (mm)	Height (mm)	Qty	Price
FSBFB501125	125	66	108	1	\$3.00
FSBFB501250	250	77	130	1	\$3.50
FSBFB501500	500	97	172	1	\$4.00
FSBFB5011000	1000	122	216	1	\$7.00
FSBFB5012000	2000	160	265	1	\$13.00



Find out more at [thermofisher.com.au](http://thermofisher.com.au)

Price valid until 31st December, 2017. © 2017 Thermo Fisher Scientific Inc. All rights reserved. Trademarks used are owned as indicated on [thermofisher.com.au/trademarks](http://thermofisher.com.au/trademarks). 149300365

**Order Placement:**

For customer service, call 1300-735-292  
To fax an order, use 1800-067-639  
To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

**Contact Us:**

For customer service, call 1300-735-292  
For service and calibration, call 1300-736-767  
To order online: [thermofisher.com.au](http://thermofisher.com.au)

**ThermoFisher**  
SCIENTIFIC