

























stomacher® bags

stomacher® 80

product code	bag type		pack quantity	pack style	paddle style	sample size	Stomacher® machine
BA6040	Standard 105mm x 155mm		2000 in sachets of 50			5-80ml	80 and 80 microBiomaster
BA6040/CLR	Closure 105mm x 155mm		200 in sachets of 10			5-80ml	80 and 80 microBiomaster
BA6040/STR	Strainer 105mm x 155mm		200 in sachets of 10			5-80ml	80 and 80 microBiomaster
BA6040/STR/DBL	Strainer double 135mm x 180mm		200 in sachets of 10			5-80ml	80 and 80 microBiomaster
BA6040/CLR/STR	Strainer closure 105mm x 155mm		200 in sachets of 10			5-80ml	80 and 80 microBiomaster
BA6040/MICRO	Micro 60mm x 185mm		50 bags			250ul-80ml	80 and 80 microBiomaster

stomacher® 3500

BA6042	Standard 380mm x 510mm		250 in sachets of 50			400-3500ml	3500 Jumbo
BA6042/HD	Heavy Duty 380mm x 510mm		200 in sachets of 10			400-3500ml (80um)	3500 Jumbo













All Stomacher® bags are made from virgin LDPE granule extruded to 65um film thickness
Irradiated sterile at 25-50kGy (Certificates available on request)

stomacher® is a registered trade mark of Seward Limited, UK













stomacher® bags

stomacher® 400 stomacher® 400

Circulator

product code	bag type		pack quantity	pack style	padding style	sample size	Stomacher® machine
BA6141	Standard 177mm x 305mm		500 in sachets of 50 BA6141/5 sachets of 5			80-400ml	Stomacher® 400 Circulator
BA6141/CPG	Standard 177mm x 305mm		500 in sachets of 50			80-400ml	Stomacher® 400 Circulator
BA6141/CLR	Closure 177mm x 305mm		250 in sachets of 10			80-400ml	Stomacher® 400 Circulator
BA6141/STR	Strainer 177mm x 350mm		200 in sachets of 10			80-400ml	Stomacher® 400 Circulator

Classic

BA6041	Standard 177mm x 305mm		500 in sachets of 50 BA6041/5 sachets of 5			80-400ml	Stomacher® 400 Classic
BA6041/CPG	Standard 177mm x 305mm		500 in sachets of 50			80-400ml	Stomacher® 400 Classic
BA6041/CLR	Closure 177mm x 305mm		250 in sachets of 10			80-400ml	Stomacher® 400 Classic
BA6041/STR	Strainer 177mm x 305mm		200 in sachets of 10			80-400ml	Stomacher® 400 Classic

Simplettes and Bag Clips

BA6051	Simplette	2000 in sachets of 50	170mm length
BA6052	Simplette	1250 in sachets of 50	250mm length
BA6099	Bag Clip	250 in sachets of 50	400ml size

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
Visit us online: www.thermofisher.com.au

©2011 ThermoFisher Scientific Inc. All rights reserved. 11-11-0100-001